

Eurofins CPT GmbH · Am Neuländer Gewerbepark 4 · D-21079 Hamburg

SATIMED UAB
attn. Ms. Evelina Marcevičienė
Sauletekio al. 7C
10257 Vilnius
LITAUEN

ProductTesting-HH@eurofins.com
www.product-testing.eurofins.com

Person in charge Ms A. Kordan - 6872
Client support Ms A. Kordan - 6872

Report date 23.10.2023
Page 1/2

Analytical report AR-23-JR-022828-01



Sample Code 799-2023-00023219

Reference	Firodenta Deep Relief Gel Forte
Client sample code	N/A
Purchase order code	N/A
Lot-no.	02/2023-1DRG
Number of received Samples	2
Ordered by	Ms. Evelina Marcevičienė
Submitted by	Ms. Evelina Marcevičienė
Carrier	UPS
Reception date	12.10.2023
Start/end of analyses	13.10.2023 / 23.10.2023

TEST RESULTS

Microbiological Analysis

LV16V	Detection of Pseudomonas aeruginosa - ISO 22717:2015 + PS/MIC/084 rev.3 2018	
	Method: EN ISO 22717, D-Cultural technique (non-chromogenic media)	
	Subcontracted to a Eurofins laboratory	
	Pseudomonas aeruginosa	Not Detected /g
LV16W	Detection of Staphylococcus aureus - ISO 22718:2015 + PS/MIC/084 rev.3 2018	
	Method: EN ISO 22718, D-Cultural technique (non-chromogenic media)	
	Subcontracted to a Eurofins laboratory	
	Staphylococcus aureus	Not Detected /g
LV16X	Detection of Candida albicans - ISO 18416:2015 + PS/MIC/084 rev.3 2018	
	Method: EN ISO 18416, D-Cultural technique (non-chromogenic media)	
	Subcontracted to a Eurofins laboratory	
	Candida albicans	Not Detected /g
LV16Y	Detection of Escherichia coli - ISO 21150:2015 + PS/MIC/084 rev.3 2018	
	Method: EN ISO 21150, D-Cultural technique (non-chromogenic media)	
	Subcontracted to a Eurofins laboratory	
	Escherichia Coli	Not Detected /g
LV16T	Enumeration and detection of aerobic mesophilic bacteria - ISO 21149:2017	
	Method: EN ISO 21149, E-Cultural technique (non-chromogenic media)	
	Subcontracted to a Eurofins laboratory	

The test results refer exclusively to the test sample provided by the customer and the scope of the tests performed.
The information about "Reference", "Client sample code", "Purchase order code", "Lot-no.", "Ordered by" and "Submitted by" were provided by the customer and may have an influence on the validity of the test results and the assessment of the results.
If a conformity statement is made, the expanded measurement uncertainty (k=2) is deducted by default when a limit value is exceeded.
Any publication of this report requires written permission. An excerpt publication is not allowed.
Eurofins CPT GmbH · Am Neuländer Gewerbepark 4 · D-21079 Hamburg
Place of execution and place of jurisdiction is Hamburg Registered Office: Hamburg - lower district court Hamburg HRB 103427 Commercial Register: Amtsgericht Hamburg HRB 103427
General Manager: Dr. Peter Schluessche
Our General Terms & Conditions of Sales are applicable - All contracts will be carried out in accordance with our General Terms and Conditions (GTC).
VAT No.: DE258239846
Bank Name: UniCredit Bank AG
BLZ: 207 300 17, Kto.-Nr.: 7000 0016 50
IBAN: DE86 2073 0017 7000 0016 50
SWIFT: HYVEDEMM17

Our General Terms & Conditions, available upon request and online at <http://www.eurofins.de/lebensmittel/kontakt/avb.aspx>, shall apply.

TAMC		<10	cfu/g
LV16U	Enumeration of yeast and mould		
	Method: EN ISO 16212, E-Cultural technique (non-chromogenic media)		
	Subcontracted to a Eurofins laboratory		
TYMC		<10	cfu/g

Signature



 Group Leader Customer Service Cosmetic/Detergent/Hygiene Products (Adelina Kordan)

The test results refer exclusively to the test sample provided by the customer and the scope of the tests performed.
 The information about "Reference", "Client sample code", "Purchase order code", "Lot-no.", "Ordered by" and "Submitted by" were provided by the customer and may have an influence on the validity of the test results and the assessment of the results.
 If a conformity statement is made, the expanded measurement uncertainty (k=2) is deducted by default when a limit value is exceeded.
 Any publication of this report requires written permission. An excerpt publication is not allowed.
 Eurofins CPT GmbH - Am Neuländer Gewerbepark 4 - D-21079 Hamburg
 Place of execution and place of jurisdiction is Hamburg Registered Office: Hamburg - lower district court Hamburg HRB 103427 Commercial Register: Amtsgericht Hamburg HRB 103427
 General Manager: Dr. Peter Schluesselche
 Our General Terms & Conditions of Sales are applicable - All contracts will be carried out in accordance with our General Terms and Conditions (GTC).
 VAT No.: DE258239846
 Bank Name: UniCredit Bank AG
 BLZ: 207 300 17, Kto.-Nr.: 7000 0016 50
 IBAN: DE86 2073 0017 7000 0016 50
 SWIFT: HYVEDEMM17

Our General Terms & Conditions, available upon request and online at <http://www.eurofins.de/lebensmittel/kontakt/avb.aspx>, shall apply.